

Reusable Temperature Probes¹

Compatible with GE Healthcare Monitors B105P/ B125P/B105M/ B125M and B155M

GE Healthcare's reusable temperature probes are intended for use with adult patients who require continuous temperature monitoring.

Features

- Provides consistent and accurate temperatures
- Dependable performance in all care areas
- Economically priced

**Our equipment.
Our accessories.
Simply compatible.**

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step



5831404

5831405

¹GE Healthcare reusable temperature probes are not available in all countries and not 510(k) cleared.

PRODUCT SPECIFICATIONS

Product use

The temperature probe is intended for use with adult patients who require continuous temperature monitoring.

Materials

Connector Plug:	Brass, Nickel Plating
Connector Premold:	PP
Connector Overmold:	TPU
Raw Cable:	TPU Jacket, Copper Conductor
Thermistor:	Copper Conductor, NTC Resistor Chip

Materials composition

Not made with any BPA, DEHP, natural rubber latex, PVC or other phthalates.

Operating environment

Temperature:	0° to 50°C (41° to 122°F)
Humidity:	10 to 95%
Pressure:	86 to 106Kpa

Storage

Temperature:	-40° to 70°C (-40° to 158°F)
--------------	------------------------------

Sterility

Temperature probes are supplied as non-sterile.

Product may not be available in all countries and regions. Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com.

Shelf Life

Shelf life is one year from date of manufacture.

Disposal considerations

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

Packaging materials

Primary: Polybag – Polyethylene (PE)

Regulatory standards

Temperature Probe Meets: ISO 10993-1, ISO 10993-5 and ISO 10993-10

Compliance

ROHS compliant

PART LIST

Part Number	Description	Quantity
5831404	Reusable Temperature Probe, Skin, Adult	1
5831405	Reusable Temperature Probe, Esophageal/Rectal, Adult	1

Distributed by:

GE Healthcare

Manufactured by:

Shenzhen Launch Electrical Co., LTD.

No. F Building

Zhonggangxing Industrial Estate

Zhangge Community, Guanlan Street

Longhua New District, Shenzhen, China



© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, the GE Monogram and Bilisoft are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

DOC2442250